



FLYING WITH YOUR

# Tandem Insulin Pump



If you plan to travel by air with your insulin pump, some advance planning may make your security screening go smoother.

Your device is safe for use during air travel and complies with CAA wireless transmission standards. It is also designed to withstand common electromagnetic interference and can be safely carried through metal detectors.

However, your Tandem Diabetes Care insulin pump should NOT be put through machines that use X-rays, including airline luggage X-ray machines and full-body scanners.

We recommend disconnecting at the infusion site and asking the security agent for an alternative screening method. If you prefer to stay connected, you can notify the agent about your pump and request to go through a standard metal detector wearing your pump.

It may be helpful to provide a printout of the information below to the security agent during the screening (this is available on [amsldiabetes.com.au/resources](https://amsldiabetes.com.au/resources) for printing).

If you have any questions about travelling with your insulin pump, please contact AMSL Diabetes Technical Support team on **1300 851 056** or at [diabetes@amsl.com.au](mailto:diabetes@amsl.com.au). We are here for you 24 hours a day, 7 days a week.

Fold Here

Show this document to the airport representative

## Dear airport representative,

This is an insulin pump from Tandem Diabetes Care, a life-sustaining medical device prescribed by a physician. It is designed to withstand common electromagnetic interference and can be safely carried through metal detectors, but it should not be exposed to any form of X-rays. This includes airline luggage X-ray machines and full-body scanners.

Please provide a screening method other than X-ray for this insulin pump. If you have any questions about this device, please call [AMSL Diabetes Technical Support on 1300 851 056](tel:1300851056)

Sincerely,  
Tandem Diabetes Care

For more information on the t:slim X2 insulin pump, please contact us on **1300 851 056** or at [diabetes@amsl.com.au](mailto:diabetes@amsl.com.au)

[amsldiabetes.com.au](https://amsldiabetes.com.au)   

© 2023 Tandem Diabetes Care, Inc. All rights reserved. ML-1000524\_D. AMSL is a Dexcom company. PR-100-267 May 2023



**Metal detectors**



**X-rays**, including full body scanners and luggage screening machines



a Dexcom company