


## The diasend<sup>®</sup> solution Interpreting the charts & graphs



*diasend<sup>®</sup> is a universal cloud-based diabetes data management system, that enables you to upload data from meters, pumps and CGMs.*

# About the diasend<sup>®</sup> solution

## diasend<sup>®</sup> solution

With the diasend<sup>®</sup> solution, data will be consolidated and viewable at [www.diasend.com](http://www.diasend.com) - without requiring any software installation. The data is presented in a clear and structured way on a secure website through graphs, tables and statistics.




## We protect your privacy


In diasend<sup>®</sup> all your private data is protected from any unauthorised access or use. All data is stored on a secure, industry standard server within the EEA region. diasend<sup>®</sup> is fully compliant with GDPR and with national privacy regulations where available. The control of sensitive data is yours and our job is to keep your data safe.

## For Clinic


1. Plug Glooko Transmitter into your power outlet



2. Connect the device to Glooko Transmitter



3. Log in at [diasend.com](http://diasend.com)



# Contents

## Glucose

<b>Glucose Standard day</b> .....	7
Patterns of high and low glucose readings over a customised period	
<b>Glucose Trend</b> .....	8
Glucose readings plotted in a graph	
<b>Glucose Meter settings</b> .....	9
<b>Before and after meal settings</b> .....	10
How to set meal times for individual patients	

## CGM

<b>CGM Standard day</b> .....	11
CGM data displayed in a box plot graph	
<b>CGM Statistics</b> .....	12
Detailed CGM statistics over a 24 hour period	

## Insulin

<b>Insulin Week</b> .....	13
Shows daily and weekly basal and bolus distribution	
<b>Insulin Bolus doses</b> .....	15
Daily bolus doses displayed in a graph	
<b>Insulin Pump settings</b> .....	16
<b>Compare pump settings</b> .....	17
View and compare pump settings from the latest uploads	
<b>Insulin Bolus adherence</b> .....	18
Log of when pump recommended boluses have been overridden	

## Comparison

<b>Comparison Logbook/table</b> .....	19
Consolidated data from insulin pumps, glucose meters and CGMs in a logbook view	
<b>Comparison Day by day</b> .....	20
Consolidated data from insulin pumps, glucose meters and CGMs displayed in tables and graphs. You can also view a PDF report with two weeks of data on one page	
<b>Comparison Day by day overview (print to pdf)</b> .....	21
Create a pdf with a 14 day overview to view on paper	

## Compilation

<b>Compilation report</b> .....	22
A summary of aggregated data from glucose meters, insulin pump and CGMs	

## PDF Wizard

<b>PDF Wizard</b> .....	23
How to set up a favorite PDF report profile	

At diasend.com data is displayed in graphs, tables and statistics. The following pages will display some of these reports and graphs you can view in diasend\*. Please note that this is a general overview of the available reports in diasend\*. All features and type of devices are not available in all countries and may therefore affect what you can view at diasend.com.


Please visit [www.diasend.com](http://www.diasend.com) for current updates regarding features and compatible devices in your specific country.

*You can find a selection of our reports in this reference guide.*

*In addition to these you can find the following reports in our system:*

- *Glucose Logbook/table*
- *Glucose Day by day*
- *Glucose Meter alarms*
- *CGM Trend*
- *CGM Day by day*
- *Insulin Trend*
- *Insulin Day by day*
- *Insulin Pump alarms*

# Better control means better healthcare



## How it works

- Both patients and healthcare professionals (HCPs) see the same reports
- Over 200 compatible devices
- Customise and generate PDF reports (print and/or save)
- Patients upload before or at Clinic
- Combine data from several devices into one patient profile

# Patients can share data with healthcare professionals (HCPs)

In addition to device data, patients can share data from activity trackers.



Share data | Export data | Subscription | **Browse & connect apps**

**IMPORTANT:** It may take a couple of hours or until next daybreak before data may appear in your diasend® reports after you have established connection. Depending on activity tracker, your activity tracker's data is then automatically pushed to diasend® at least once a day when new data is available.

**Fitbit**  
Fitbit offers wearable devices to help you lead a healthier, more active life.  
[Connect](#)

**Runkeeper**  
The Runkeeper app is a personal trainer in your pocket.  
[Connect](#)

**Dexcom (GS Mobile / G6 app)**  
The Dexcom system uses a sensor and app to continuously track your glucose levels and trends for better diabetes management.  
[Connect](#)

**Abbott FreeStyle LibreLink**  
*Available in Sweden and Norway only.*  
[Read more](#)  
The FreeStyle Libre sensor is read via the FreeStyle LibreLink app on your Android smartphone or iPhone. Connecting your LibreView account to your diasend® account means that your data will continuously sync with diasend®.  
[Connect](#)

**Medtronic**  
A6 TouchCare® combines a Patchpump and a CGM to communicate to the PDM and then to the EasyTouch app (or EasySense app if you use the CGM standalone) on iOS or Android smartphones. When you connect your EasyTouch/EasySense account to diasend®, the data will be continuously and automatically uploaded to your account.  
[Connect](#)

**mylife™ App / mylife™ Cloud**  
Log in with your mylife Cloud credentials and connect it to your diasend® account. Have your therapy data synchronized at all times, it does not matter whether you have recorded your data via the mylife App or the mylife Software.  
[Connect](#)

**Eversense**  
The Eversense Mobile Application works with the Eversense Continuous Glucose Monitoring System. The mobile app provides easy access to CGM data and eliminates the need to carry a dedicated receiver device.  
[Connect](#)

*Patients can easily connect their favorite activity tracker to their diasend® account.*

Patients can register for free, upload from home and share their data with their healthcare professionals at [www.diasend.com/register](http://www.diasend.com/register)

**Physical activity targets**  
Set your target values

Steps per day	Calories (kcal) to burn per day
<input type="text" value="10000"/>	<input type="text" value="2500"/>

Activity	
Avg steps / day	Avg kcal / day
<b>4047</b> steps	<b>950</b> kcal
40% of 10000 (target)	38% of 2500 (target)

# Report information

*diasend<sup>®</sup> gives you increased accessibility to glucose readings, insulin doses and CGM data. This gives the user the ability to become more involved in their diabetes management and care. diasend<sup>®</sup> is easily customised, giving the individual patient and HCP only the data they are looking for. Our customisation ultimately allows for a more effective conversation between the patient and the HCP.*

## **Compilation Report**

Allows the user to get an overall picture of the data uploaded for that time period.

## **Glucose Standard Day Report**

Can be used to identify patients who are not testing regularly and/or look at different time periods to identify problems.

## **CGM Reports**

Clearly display data to help assist the HCP in identifying problem areas and/or times for each patient. You can use the CGM Standard Day report to identify the range of CGM readings by the time of day or look at the CGM statistics to check standard deviation and averages by the hour.

## **Day by Day Report**

Easily compare different hours of the day, or weekdays against weekends. This may help the patients identify activities or events that have influenced their values.

## **Insulin Bolus Dose Report**

Check for patterns of how the patient is bolusing. This will help identify when the patient is doing well and where they need to improve.

## **Bolus Adherence Report**

Easily see what the pump has calculated for the patient's bolus against what the patient actually delivered.

## **Insulin Pump Settings Report**

View historical data to easily compare and contrast pump settings from different time periods.

## **The Comparison Logbook/table Report**

Gives the HCP insight into glucose measurement, carb intake and cannula fills, as well as boluses, priming, and suspend events, which assists in evaluating events by time and day (and may aid in identifying patterns). This report view also displays ketones if this is saved in the meter.

## **Settings Tab**

Easily change the default blood glucose (BG) target range for the clinic and also individual patients as well as setting up customised PDF reports.


# Glucose

## STANDARD DAY

# Glucose standard day

An overview of glucose readings plotted by time during a standard 24 hour day.

Can be used to identify patients who are not testing regularly and/or look at different time periods to identify problems.




- 1 You can select a midnight-to-midnight or noon-to-noon view over the selected date range.
- 2 You can choose to show lines to connect the values relating to each day of the week in different colors.
- 3 View the mean value by clicking on this link.
- 4 Statistics can be found at the bottom of each report.

# Glucose trend

This displays a trend overview of glucose readings by date.

1 This view also offers the possibility to look at values during specific time periods such as before or after meals.

2 Click here to Show lines.



# Glucose meter settings

## Glucose

### METER SETTINGS

Glucose CGM Insulin Comparison Compilation My profile

LOGBOOK/TABLE STANDARD DAY TREND DAY BY DAY METER ALARMS **METER SETTINGS** A1C

Meter settings for serial number: 1 Print to PDF

General		Correction target (Range)		Correction factor		I:C (g option)	
Setting	Value	Setting	Value	Setting	Value	Setting	Value
Food unit	Grams of carbs	All day	5 mmol/L - 8 mmol/L	All day	2.7 mmol/L	Morning	10 g
Dose Increments	1 U					Midday	12 g
Correction Target Type	Range					Evening	14 g
Insulin Log Feature	Off						
Insulin Calculator	Advanced						
Notes Feature	Off						

General		Health events	
Setting	Value	Setting	Value
Meal excursion	4.2 mmol/L	Exercise 1	-25 %
Snack limit	20 g	Exercise 2	-42 %
Active timeout	02:45:00	Stress	22 %
Offset timeout	01:30:00	Illness	33 %
		Premenstrual	16 %

Timeblock					
Setting	00:00-06:30	06:30-12:30	12:30-18:30	18:30-22:30	22:30-00:00
Target interval min	3.8 mmol/L	4 mmol/L	4.4 mmol/L	3.5 mmol/L	4.4 mmol/L
Target interval max	6.6 mmol/L	8.1 mmol/L	8.2 mmol/L	8.5 mmol/L	6.7 mmol/L
Carb ratio, insulin	1.9 U	2.6 U	2.1 U	2.5 U	2 U
Carb ratio, carbs	11 g	13 g	14 g	15 g	12 g
Insulin sensitivity, insulin	0.9 U	0.7 U	3.1 U	1.8 U	50 U
Insulin sensitivity, carbs	2.9 mmol/L	0.4 mmol/L	54.3 mmol/L	3.3 mmol/L	36 mmol/L

This shows the settings for meters that have advanced settings activated.

Meter settings data is displayed in different formats depending on the device that is uploaded.

Examples of meter settings are shown on the left.

# Before and after meal settings

You can choose if you want to display time intervals, or you can set it to show before and after meal times. You will see this in the Compilation (page 21) and Glucose trend (page 8) reports.

- 1 Click on the patient profile tab.
- 2 Expand the registration form to set the meal times for the patient.
- 3 Click on permanent intervals.
- 4 Set the meal times.
- 5 Click to save the information.

Glucose CGM Insulin Comparison Compilation

(\*) means that the fields are compulsory

Expand registration form Print to PDF

My profile

Personal details

First name (\*) Last name (\*) Personal ID Date of birth (yyyy-mm-dd) Gender

Jane Doe [ ] 1979 Dec 20 Female

Intervals

There are no permanent intervals for the patient (default setting)  
Explanation: this option means that there are no specific intervals, such as before breakfast.

Permanent intervals  
Explanation: this option means that settings for specific intervals will be made, such as before breakfast.  
Please note: Placing intervals close to each other than the selected interval will cause problems.

The intervals durations before meals: 1  
The intervals durations after meals: 1

	Weekdays	Weekends
When do you have breakfast?	00	00
When do you have lunch?	00	00
When do you have dinner?	00	00
When do you go to bed?	00	00


Save my information

When you're done, click "Save my information" below.

Save my information

# CGM standard day


## CGM STANDARD DAY




*This mode shows the distribution of values per hour of the day, spanning over multiple days.*

- 1** The median line represents the middle of the value series within each hour.
- 2** The darker blue area represents the Interquartile range ("IQR"), meaning the middle 50% of the values are within this range.
- 3** The lighter blue area represents the 10th and 90th percentile, meaning that middle 80% of the values are within this range.
- 4** The maximum and minimum values are also visualised.
- 5** You can click on this link to choose to see the values as an AGP, Box plot or modal graph.

Box plot



Modal



# CGM statistics

This table includes detailed CGM statistics over the time period selected by hour of the day, check SD (standard deviation) and averages by the hour.

- 1 Statistics from a particular hour of the day.

**CGM**                      **STATISTICS**

Glucose   **CGM**   Insulin   Comparison   Compilation   My profile

STANDARD DAY   TREND   DAY BY DAY   **STATISTICS**

Period: 08/19/2014 - 09/01/2014, 14 days    ▼

Include manually entered records

▼

	00	01	02	03	04	05	06	07	08	09	10	11
# of CGM Readings	159	167	168	168	168	168	168	168	168	162	156	152
Median CGM Value (mmol/L)	8.3	8.9	9.4	11.2	12.6	12.8	13.8	12.4	9.6	7.9	7.2	6.7
Avg CGM Value (mmol/L)	8.7	8.9	9.9	10.9	12.6	13.5	13.4	12.8	11.5	9.8	8.4	7.3
Min (mmol/L)	3.4	3.2	4.8	4.3	5.2	5.1	4.3	3.7	3.3	3.9	3.6	3.4
25% Quartile	6.2	7.1	7.5	8.4	10.1	10.7	10.4	9.6	7.3	6.2	5.8	5.1
75% Quartile	10.9	10.4	11.4	12	13.9	15.7	14.9	16.1	16.9	12.8	9.9	8.7
Max (mmol/L)	15.5	14.5	17.1	20.5	20.7	22.3	22.3	22.3	22.3	22.3	20.2	15.9
SD (mmol/L)	2.8	2.6	2.9	3.4	3.5	4.1	4.7	5	5.5	5.1	3.7	2.7

	12	13	14	15	16	17	18	19	20	21	22	23	Totals
# of CGM Readings	137	129	132	132	143	139	145	156	155	143	132	143	3658
Median CGM Value (mmol/L)	7.2	8.3	9.9	9.4	9	7.6	6.6	8.6	9.8	10	8.6	7.5	9.4
Avg CGM Value (mmol/L)	7.8	9.1	10.4	10.3	10.4	8.9	7.1	9.1	10.3	10.7	9.6	8.4	10.1
Min (mmol/L)	3.3	5.2	6.8	6.6	5.6	4.9	2.8	4.3	5.1	4.4	4.2	3.7	2.8
25% Quartile	6.3	7.2	8.7	8.5	7.8	6.1	4.9	7.1	7.6	7.7	6.1	5.2	7.1
75% Quartile	10.2	10.5	12.1	12.4	12.6	10.8	7.7	10.9	13.2	12.9	11	10.4	12.4
Max (mmol/L)	11.6	13.9	15.7	19.1	20.8	20.2	15.1	14.7	17.2	19.3	19.2	19.1	22.3

# Insulin week




Daily and weekly basal and bolus distribution.

- 1** Daily basal-bolus numbers can be viewed in the bar graph.
- 2** You can view the insulin distribution for the week in the pie chart.



# Insulin bolus doses



*Bolus doses displayed in graph format to easily view daily bolus activity by time of day.*

*This graph gives you a general overview of the bolus doses over a period of time. Check for patterns of how the patient is bolusing. This will help identify when the patient is doing well and where they need to improve.*

*For example you can see gaps of missed boluses.*

- 1** A combo bolus may be used for high fat/carb meals providing a percentage of the dose immediately and then a slow infusion of insulin spread out over a set amount of time.

# Insulin pump settings

## Insulin

## PUMP SETTINGS

This provides the information of the current and historical settings in the pump to easily review, compare and adjust as necessary.

- 1 Select and view the pump settings from every upload.
- 2 Print a comparison of the latest available pump settings (see page 17).
- 3 Select one of these options to print the displayed pump setting in an expanded version, or in a minimised version where all pump settings are compressed onto 1 page (see page 22 for an example).
- 4 View the Bolus, Basal, General and CGM settings. In this report you can also view I:C ratio, ISF and glucose target ranges.

The screenshot shows the 'Insulin' tab selected in the top navigation bar. Below it, the 'PUMP SETTINGS' sub-tab is active. A period selector shows '08/19/2014 - 09/01/2014, 14 days'. A dropdown menu for 'Insulin pump settings for serial number:' is open, showing '09/01/2014 08:01 (US/Eastern)'. Three print buttons are visible: 'Print comparison of pump settings to PDF', 'Print to PDF on 1 page', and 'Print to PDF'. The settings are organized into three columns: Bolus, General, and CGM Settings. The 'Basal' section is highlighted with a red circle.

Setting	Value
Audio Bolus Enable	Disabled
Audio Bolus Stepsize per program keypress	0.5 U
Advanced Bolus Options enable	Enabled
Bolus Reminder Options enable	Enabled
Bolus Delivery Speed	Normal
Bolus	20 U

Setting	Value
Language Selection Index	English
Last Keypress to display timeout	30 s
Auto-Off Enable	Enabled
Auto-Off Timeout	16 h
Max 2-Hr limit	30 U
Occlusion Sensitivity Level	High
Insulin-On-Board	Enabled
Insulin-On-Board Duration	4 h
Sick days, BG over limit	17 mmol/L
Sick days, check ketones	2 h
Sick days, check BG	2 h
Low Cartridge Warning Level	40 U

Setting	Value
Transmitter Sound Level	Vibration
Other Sound Level	Vibration
Glucose High Alert Limit	7.8 mmol/L
Glucose Low Alert Limit	3.9 mmol/L
Glucose Rise Alert Limit	0.1 mmol/L
Glucose Fall Alert Limit	0.1 mmol/L
Glucose Low Alert Snooze Time	60 min
Glucose High Alert Snooze Time	30 min
Transmitter Out of Range Alert Snooze Time	60 min
Glucose Low Enable	Enabled
Glucose High Enable	Enabled
Glucose Rise Enable	Enabled

Setting	Value
Max Basal	2 U/h
Max Total Daily Dose	80 U
Active basal program	1

# Compare pump settings

(PRINTED COMPARISON OF PUMP SETTINGS)

Patient: Jane Doe  
 Patient ID:   
 Date interval: 08/19/2014 - 09/01/2014  
 Number of days: 14  
 Glucose meters:   
 Insulin pump:   


## Insulin: Pump settings comparison

Insulin pump settings for serial number: XX-XXXX-XX (Animas Vibe). The report compares the settings from the latest eight uploads. Changes are marked.

	Upload date: 01/09/2014	03/05/2014	08/04/2014	08/06/2014	08/08/2014	08/18/2014	08/19/2014	09/01/2014
<b>General</b>								
Language Selection Index	English	English	English	English	English	English	English	English
Last Keypress to display timeout	30 s	30 s	30 s	30 s	30 s	30 s	30 s	30 s
Auto-Off Enable	Enabled	Enabled	Enabled	Enabled	Enabled	Enabled	Enabled	Enabled
Auto-Off Timeout	16 h	16 h	16 h	16 h	16 h	16 h	16 h	16 h
Max 2-Hr limit	30 U	30 U	30 U	30 U	30 U	30 U	30 U	30 U
Occlusion Sensitivity Level	High	High	High	High	High	High	High	High
Insulin-On-Board	Enabled	Enabled	Enabled	Enabled	Enabled	Enabled	Enabled	Enabled
Insulin-On-Board Duration	4 h	4 h	4 h	4 h	4 h	4 h	4 h	4 h
Sick days, BG over limit	17 mmol/L	17 mmol/L	17 mmol/L	17 mmol/L	17 mmol/L	17 mmol/L	17 mmol/L	17 mmol/L
Sick days, check ketones	2 h	2 h	2 h	2 h	2 h	2 h	2 h	2 h
Sick days, check BG	2 h	2 h	2 h	2 h	2 h	2 h	2 h	2 h
Low Cartridge Warning Level	40 U	40 U	40 U	40 U	40 U	40 U	40 U	40 U
Time format	24 h	24 h	24 h	24 h	24 h	24 h	24 h	24 h
BG unit	mmol/L	mmol/L	mmol/L	mmol/L	mmol/L	mmol/L	mmol/L	mmol/L
<b>Bolus</b>								
Audio Bolus Enable	Disabled	Disabled	Disabled	Disabled	Disabled	Disabled	Disabled	Disabled
Audio Bolus Stepsize per program keypress	1.0 U	0.5 U	0.5 U	0.5 U	0.5 U	0.5 U	0.5 U	0.5 U
Advanced Bolus Options enable	Enabled	Enabled	Enabled	Enabled	Enabled	Enabled	Enabled	Enabled
Bolus Reminder Options enable	Enabled	Enabled	Enabled	Enabled	Enabled	Enabled	Enabled	Enabled
Bolus Delivery Speed	Normal	Slow	Normal	Normal	Normal	Normal	Normal	Normal
Max Bolus	20 U	20 U	20 U	20 U	20 U	20 U	20 U	20 U

- 1 The date of each upload.
- 2 For ease of reference all changes, compared to the latest uploads to diasend®, are highlighted in the report.

# Insulin bolus adherence

Log of when bolus overrides have occurred. This report allows you to see what the pump calculated as the patient's bolus against what the patient actually delivered.

The report displays:

- 1 Insulin actually delivered.
- 2 Amount of insulin suggested by bolus calculator.
- 3 Bolus Type.
- 4 Pre-bolus BG reading.
- 5 Post-bolus BG reading

**Insulin**

**BOLUS ADHERENCE**

Glucose CGM **Insulin** Comparison Completion My profile

WEEK TREND DAY BY DAY BOLUS DOSES **BOLUS ADHERENCE** PUMP ALARMS PUMP SETTINGS

Period: 08/19/2014 - 09/01/2014, 14 days Select time interval ▼

Print to PDF ▼


14 deviating boluses detected out of 100 bolus calculator-assisted boluses (out of 100 boluses in total)

Date	Time	Delivered (U)	Calculated (U)	Bolus Type	Duration (min)	Pre-Bolus BG (mmol/L)	Post-Bolus BG (mmol/L)
08/19/2014	14:07	2.00	0.00	ezBG		08/19/2014 14:07 : 13.0	08/19/2014 15:33 : 10.2
08/21/2014	11:51	4.75	5.75	ezCarb (Normal)		08/21/2014 11:51 : 3.8	08/21/2014 11:56 : 4.0
08/21/2014	11:56	6.90	0.90	ezCarb (Normal)		08/21/2014 11:56 : 4.0	08/21/2014 14:36 : 6.8
08/22/2014	02:45	1.55	1.80	ezBG		08/22/2014 02:45 : 12.1	08/22/2014 03:42 : 12.5
08/22/2014	03:42	1.20	1.50	ezBG		08/22/2014 03:42 : 12.5	08/22/2014 07:42 : 8.3
08/22/2014	07:42	0.85	0.80	ezBG		08/22/2014 07:42 : 8.3	08/22/2014 15:33 : 9.9
08/22/2014	15:33	1.00	0.25	ezBG		08/22/2014 15:33 : 9.9	08/22/2014 19:54 : 7.3
08/22/2014	19:54	0.40	0.00	ezBG		08/22/2014 19:54 : 7.3	08/23/2014 03:43 : 4.4
08/23/2014	03:43	1.10	0.00	ezBG		08/23/2014 03:43 : 4.4	08/23/2014 07:21 : 10.2
08/23/2014	12:31	2.05	0.00	ezBG		08/23/2014 12:31 : 10.6	08/23/2014 15:54 : 18.9
08/23/2014	15:54	6.25	5.90	ezBG		08/23/2014 15:54 : 18.9	08/23/2014 16:41 : 20.5
08/26/2014	22:20	4.00	0.00	ezBG		08/26/2014 22:20 : 18.4	08/27/2014 13:14 : 9.8
08/31/2014	08:39	7.00	6.75	ezBG		08/31/2014 08:39 : 20.0	08/31/2014 09:34 : 22.2
09/01/2014	12:37	7.90	7.00	ezCarb (Normal)		09/01/2014 12:37 : 7.0	

# Comparison logbook/table

## Comparison

### LOGBOOK/TABLE



Consolidated data from insulin pump and glucose meters in a logbook view. Gives the HCP insight into glucose measurement, carb intake and cannula fills, as well as boluses, priming, and suspend events, which assists in evaluating events by time and day (and may aid in identifying patterns). This report view also displays ketones if this is saved in the meter.

- 1 The free standing numbers are glucose readings.
- 2 "g" stands for grams of carbohydrates.
- 3 "U" stands for units of insulin.
- 4 The box indicates that a correction was made.
- 5 You can view pump events.
- 6 Red = above target range  
Blue = below target range  
Green = within target range
- 7 Daily averages also available


# Comparison day by day

## Comparison

### DAY BY DAY

A day by day view of consolidated data from insulin pumps, glucose meters and CGMs in table and graphs.

- 1 Click on this button to see a 2 week overview on one page (see page 21 for a sample).
- 2 Details of events and alarms can be seen if you hover over the icons above the graph.
- 3 The bolus and basal graph will display insulin data which includes basal rate, temporary basal rate, boluses doses, combination boluses.
- 4 You can view the daily total basal and bolus insulin distribution which is displayed in a separate pie chart.
- 5 You can view the carbs in this graph as well as the CGM curve and calibrations if CGM has been uploaded.



# Comparison day by day overview (print to pdf)

1




Day by day overview ▾

Create a pdf with a 14 day overview!


All reports in diasend<sup>®</sup> can also easily and quickly be exported and printed in PDF by setting up the PDF Wizard. Please see page 23 for information on personalisation/customisation of PDF reports.

1

By clicking here you can create a day by day overview containing several weeks. This gives you the ability to compare data from one week to the next.



## Comparison : Day by day overview



# Compilation report

A summary of aggregated data from glucose meters, insulin pumps and CGMs. Allows the user to get an overall picture of the data uploaded for that time period.

- 1 View average BG and SD by time of day.
- 2 View detailed CGM data such as average by time of day and AUC (area under the curve).
- 3 View detailed insulin pump and carb information such as average days between cannula fills and average days between primes.

## Compilation

Glucose CGM Insulin Comparison **Compilation** My profile

Period: 08/29/2014 - 09/11/2014, 14 days Select time interval


Glucose		CGM		Insulin		Carbs		Activity	
Average <b>11.2</b> mmol/L		Average <b>10.1</b> mmol/L		Average daily dose <b>29.2 U</b>		Average carbs / day <b>67 g</b>		Avg steps / day <b>4047</b> steps	Avg kcal / day <b>950</b> kcal
SD = 4.7	# = 46	SD = 4.1	# = 3658	SD = 26	# days = 24	SD = 117	# = 26	40% of 10000 (target)	38% of 2500 (target)
Avg # / day = 1.9		Avg # / day = 152.4		Avg # bolus doses/day = 4.9		Avg # / day = 1.1			

Interval	Avg BG	# BG	SD
00:00-06:00	12	4	5.2
06:00-08:00	11.1	3	2.8
08:00-10:00	16.4	6	4.4
10:00-12:00	6.6	3	3.8
12:00-14:00	9.1	10	2.1
14:00-16:00	11.4	8	3.4
16:00-18:00	13.6	2	6.9
18:00-20:00	11.2	2	3.9
20:00-22:00	9.5	7	3.5
22:00-24:00	18.4	1	0

### Glucose (mmol/L)

Glucose values summary	
Average (mmol/L)	11.2
Median (mmol/L)	10.6
Highest value (mmol/L)	22.2
Lowest value (mmol/L)	3.8
Standard deviation (SD)	4.7
Values per day	1.9
Number of values	46
Values above goal (10 mmol/L)	26
Values within goal (3.9-10 mmol/L)	19
Values below goal (3.9 mmol/L)	1


Interval	Avg BG	# BG	SD
00:00-06:00	12	4	5.2
06:00-08:00	11.1	3	2.8
08:00-10:00	16.4	6	4.4
10:00-12:00	6.6	3	3.8
12:00-14:00	9.1	10	2.1
14:00-16:00	11.4	8	3.4
16:00-18:00	13.6	2	6.9
18:00-20:00	11.2	2	3.9
20:00-22:00	9.5	7	3.5
22:00-24:00	18.4	1	0



### CGM (mmol/L)

CGM readings summary	
Average (mmol/L)	10.1
Median (mmol/L)	9.4
AUC high > 10 mmol/L	1.6
AUC low < 3.9 mmol/L	0
Highest value (mmol/L)	22.3
Lowest value (mmol/L)	2.8
Standard deviation (SD)	4.1
Values per day	152.4
Number of values	3658
Values above goal (10 mmol/L)	1599
Values within goal (3.9-10 mmol/L)	2006

Interval	Avg	#	SD
00:00-06:00	10.8	999	3.8
06:00-08:00	13.1	338	4.9
08:00-10:00	10.7	334	5.4
10:00-12:00	7.9	312	3.3
12:00-14:00	8.4	270	2.4
14:00-16:00	10.4	271	2.5
16:00-18:00	9.7	284	4
18:00-20:00	8.1	305	3
20:00-22:00	10.4	303	3.7
22:00-24:00	9	277	4.3




AUC high > 10 mmol/L	1.6
AUC low < 3.9 mmol/L	0

### Insulin

Insulin doses summary	
Average daily insulin (U)	29.2
Standard deviation (SD)	25.8
Average daily basal (U)	9.7
Average daily bolus (U)	19.5
Average bolus doses/day	4.9
Average days between cannula fills	2.3
Average days between primes	2.2

Bolus calculation summary	
Avg # ezBG Boluses/day	1.1 (23%)
Avg # ezCarb Boluses/day	0.8 (15%)
Avg # Combo Boluses/day	0.1 (2%)
Avg # Normal Boluses/day	3 (61%)
Bolus overrides/total boluses	12%
Avg # bolus overrides/day	0.6
Avg # bolus ezBG overrides/day	0.5
Avg # bolus ezCarb overrides/day	0.1
Avg # carbs/ezCarb Bolus	58 g
Avg # Insulin Units/ezCarb Bolus	6



Average days between cannula fills	2.3
Average days between primes	2.2

# PDF wizard

PDF Wizard

This wizard will allow you to create a single PDF file with all reports of your choice.

Select the reports you would like to include:

End date **1**  Use same time span for all?  No  Yes

Include manually entered records

Color **2**  Black and white **3**

**Glucose**

Logbook/table

Standard day

Trend

Day by day

Meter alarms

Meter settings

**CGM**

Standard day

Trend

Day by day

Statistics

**Insulin**

Week

Trend

Day by day

Bolus doses

Bolus adherence

Pump alarms

Pump settings

Pump settings comparison

Options

Options

Options

Two weeks

**4**

EMR PDF

PDF

★ Dr. Huffman

Julie

Go to PDF Wizard

PDF Settings

**5**

SETTINGS ACCOUNT PROFILE LOG OUT

EMR/PDF Wizard Password Default target range

The EMR/PDF Wizard allows you to create Favorite Reports Profiles, which allow you to consolidate multiple reports into a single Adobe PDF file from record systems that are compatible with PDF files.

Favorite Reports Profiles:

Name of Profile	Default	Action
	<input type="checkbox"/>	Edit  Delete
Dr. Huffman	<input checked="" type="checkbox"/> ★	Edit  Delete
Julie	<input type="checkbox"/>	Edit  Delete

[+ Add new Favorite Reports Profile](#)

Action

When I click the EMR/PDF Wizard button, I would like diasend® to perform the following action:

PDF ★ Instantly generate the default PDF profile that has been selected in the settings/pdf wizard tab

PDF Always guide me through the EMR/PDF Wizard.

PDF Always let me choose among favorite profiles already set up or guide me through the EMR/PDF wizard.

Allows customisation and management of report preferences consolidated into one PDF file. This includes the ability to add favorite profiles per user.

- 1** Set the end date of the desired time interval.
- 2** The PDF Wizard allows you to choose which reports to include.
- 3** Set the time span.
- 4** In the patient list you can click on the PDF icon to create a compiled PDF report or click on the downward pointing arrow to select the pre-set profile you wish to view/print.
- 5** You can also find the PDF Wizard under Settings in your account. There you can also add new customised report profiles.

*diasend® is indicated for use by individuals or healthcare professionals in the home or health care facilities for transmitting data from home monitoring devices such as glucose meters and insulin pumps to a server database to support diabetes management.*

For more information on diasend by Glooko please contact us on **1300 851 056** or at [diabetes@amsl.com.au](mailto:diabetes@amsl.com.au)

[amsldiabetes.com.au](http://amsldiabetes.com.au)    

To contact Glooko directly, visit [glooko.com/int](http://glooko.com/int) or email [support@diasend.com](mailto:support@diasend.com). PR-100-49 March 2020

