

## Dexcom G4® Receiver Icons

---



This icon means you need to enter a calibration. This may occur when it is time for your 12 hour calibration update or any other time an additional calibration is needed



This icon means you need to enter one more blood glucose value in order to calibrate the system and start receiving sensor glucose readings



This icon means the receiver does not understand the sensor signal temporarily. Wait for more prompts and do not enter any blood glucose values



This icon means the receiver and sensor/transmitter are not communicating. Make sure the receiver and sensor/transmitter are within 6m of each other



This icon means the sensor is not calibrating correctly. Enter at least one more calibration blood glucose value after about 10-15 minutes



This icon means the sensor is not calibrating correctly. Wait a minimum of one hour and then enter one more blood glucose value



This icon means the receiver has detected a potential problem with the sensor signal. Do not enter any blood glucose values when you see this. Wait 30 minutes for more prompts