


# t:slim X2™ Insulin Pump Infusion Sets & Cartridge




## inset™ II



- 90° insertion angle
- Integrated inserter
- Teflon cannula

## inset™ 30



- 30° insertion angle
- Integrated inserter
- Teflon cannula

## contact™ detach



- Manual straight in insertion
- Stainless steel/metal cannula

## comfort™ set



- Manual angled set
- Teflon cannula

## Infusion Set NDSS Order Numbers

Inset II Grey 6mm/60cm, 10	545
Inset II Grey 9mm/60cm, 10	546
Inset II Grey 6mm/110cm, 10	547
Inset II Grey 9mm/110cm, 10	548
Inset 30 Grey 13mm/110cm, 10	537
Inset 30 Grey 13mm/60cm, 10	540
Contact Detach 6mm/60cm, 10	673
Contact Detach 6mm/110cm, 10	674
Contact Detach 8mm/60cm, 10	675
Contact Detach 8mm/110cm, 10	676

Comfort Full Set 17mm/60cm, 10	541
Comfort Full Set 17mm/110cm, 10	543
Comfort Full Set 13mm/60cm, 10	549
Comfort Full Set 13mm/110cm, 10	551

## t:slim X2™ Cartridge

Tandem t:slim™ Cartridge Luer Lock	720
------------------------------------	-----

## Cannula Fill and Site Change Guidelines for Infusion Sets

	Fill	Site Change
<b>Inset II</b>	6mm: 0.3 units, 9mm: 0.5 units	2-3 days
<b>Inset 30</b>	0.7 units	2-3 days
<b>Contact Detach</b>	0 units/no fill	1-2 days
<b>Comfort Set</b>	0.7 units	2-3 days

For more information on the t:slim X2 Insulin Pump, please contact the AMSL Diabetes Customer Care Team on **1300 851 056**.

[amsldiabetes.com.au](http://amsldiabetes.com.au)    

Always read the label and use only as directed. Consult a healthcare professional to see which product is right for you.

ARTG 116234, ARTG 304681. PR-100-50

