

# Dexcom Clarity User Account

## Dexcom G5 Mobile System

1. Visit <https://clarity.dexcom.com/>
2. Select 'Dexcom CLARITY for Home Users'



3. The Username and Password is your Dexcom G5 Mobile app Username and Password. *You must either recover it via email or remember upon set up of Dexcom app.* Dexcom G4 PLATINUM users will need to create a new account by clicking 'Create a Dexcom Account'.
4. If you own a Receiver, you can upload it by downloading the uploader and plugging the Receiver into the computer.
5. Dexcom CLARITY receives streamed glucose data from your Dexcom G5 Mobile CGM System about 3 hours after the data is recorded.

**NOTE: Users do not need to upload their Dexcom G5 Mobile as it is constantly updated in their account once they log in.**

