


# An introduction to



## Quickly Asses Key Metrics



## Seamless Sharing

Patients using the Dexcom® continuous glucose monitoring (CGM) app automatically and continuously send their glucose data to their CLARITY® accounts.


Patients authorise data sharing using the CLARITY Clinic sharing code to link the patient and clinic accounts.

Encourage patients to use the CLARITY app so they can receive weekly notifications of progress.



## Relevant Glucose Insights

CLARITY helps prioritise clinical issues based on discussions with the patient.



Register your clinic at [clarity.dexcom.eu/professional](http://clarity.dexcom.eu/professional)


\*Estimated A1C is only an estimate and does not replace the HbA1c test. Dexcom CLARITY needs at least 12 days of sensor use with 80% or greater wear time to provide an estimated A1C. If "N/A" is displayed here, there is not enough data to estimate A1C.

# Dexcom CLARITY provides relevant patterns and trends to help with diabetes management.

Offering a full range of reports to suit your needs, including patterns, trends, and statistics.


Analyse individual days and isolated events.

## Daily




Visualise patterns using all CGM data points.

## Overlay



View the variability around the mean glucose and patterned areas of highs and lows.

## AGP (Ambulatory Glucose Profile)



Encourage progress and evaluate changes made.

## Compare

