


ONETOUCH

Glucose Meter Comparison

Always read the label and use only as directed. Read the warnings available on www.amsldiabetes.com.au/ resources before purchasing. Consult your healthcare professional to see which products are right for you.

Product Name	OneTouch Verio IQ®	OneTouch Verio Flex®
Glucose Meters	 <p>SMART.</p>	 <p>SIMPLE.</p>
Technology	PatternAlert™ Technology helps notify you of high or low glucose patterns with easy one step meal tagging*	ColourSure™ technology instantly shows when blood glucose results are in or out of range
Blood Sample	0.4 µL	0.4 µL
Test Time	5 seconds	5 seconds
Compatible Software	Glooko diasend web-based software and app	Glooko diasend web-based software and app
Test Memory	Stores up to 750 tests	Stores up to 500 tests
Battery	Re-chargeable with OneTouch mini USB	Long lasting lithium battery
Test Strip Light	Illuminated test strip port	N/A
Display	Colour screen	Black and White
Consumables	Strips: Verio Test Strips (NDSS Code: 120) Lancets: OneTouch UltraSoft Lancets	Strips: Verio Test Strips (NDSS Code: 120) Lancets: OneTouch Verio Flex Delica Lancets
Customer Service	Australian-based customer care (7am - 6pm Monday to Friday)	Australian-based customer care (7am - 6pm Monday to Friday)

For more information on the OneTouch product range or to order a free glucose meter please contact us on **1300 851 056** or at diabetes@amsl.com.au

amsldiabetes.com.au    



*The meter uses results tagged "before meal" to find high patterns. Tagging is not necessary to find low patterns.
OneTouch Verio Flex and OneTouch Verio IQ are registered trademarks of LifeScan. ARTG 168174, PR-100-194 July 2019