

## Step 1: How to create a patient account

You will receive an e-mail from your nurse/doctor with the title "Your healthcare professional invites you to join diasend<sup>®</sup>".

Click on the link in this e-mail that says [Click here to get started](#) and follow the steps below.

If you haven't received an e-mail, you can create an account at [www.diasend.com](http://www.diasend.com) by clicking the "Register here" button.

### 1 Login details

Enter your email address, which will be your username, as well as a password of your choosing for your diasend<sup>®</sup> personal account. Click continue.

### 2 Personal information

Enter your personal information then click continue.

### 3 Share data

If you received an invite from your healthcare professional via e-mail, the Clinic ID number will already be filled in. This is the unique ID number for your clinic which makes it possible for you to share your data with them. Click continue.

If you didn't receive an invite, you can add this Clinic ID later. Finish the registration, then log in to your account and go to the "Account admin" tab.

### 4 Confirm registration

Check the box to approve terms and conditions and confirm the registration. Click continue.

### 5 Download and install diasend<sup>®</sup> uploader

Install diasend<sup>®</sup> uploader on your PC or Mac, by clicking on the relevant icon. Follow the installation instructions on the screen.

Need support?

Contact AMSL Diabetes on 1300 851 056 or at [diabetes@amsl.com.au](mailto:diabetes@amsl.com.au)

## Step 2: How to start uploading data to diasend®

1

### Start diasend® Uploader


Double click on the diasend® uploader icon on your computer (if you have a Mac please go to the "applications folder" or "launchpad" to find diasend® uploader).



2

### Connect the device cable


Connect the cable for your device to a USB port on your computer.



3

### Connect your device to the cable


Connect the device to the cable or place it in front of an IR dongle.



4

### Device will start to upload


Do not unplug your device until the uploading process is completed.



5

### Verify your account


Enter your username (entire email address and the password you created when registering for your account). Click Sign in.



6

### View your data

Log in at [www.diasend.com](http://www.diasend.com) and view your data.



7

### Access your data from your mobile device

Download the diasend® Mobile App from Google Play (Android) or the App Store (iOS) to review and upload your data from your mobile device.



**This software does not provide any medical advice and should not be used for that purpose. Consult your healthcare professional for any clinical advice and to see if this product is right for you.**

For more information on diasend by Glooko please contact us on **1300 851 056** or at [diabetes@amsl.com.au](mailto:diabetes@amsl.com.au)

[amsldiabetes.com.au](http://amsldiabetes.com.au)

