

dexcomG6

Continuous Glucose Monitoring System

**MAKE
KNOWLEDGE
YOUR
SUPERPOWER**



Always read the label and use only as directed. Read the warnings available on amsldiabetes.com.au/resources before purchasing. Consult your healthcare professional to see which products are right for you.

INTRODUCING THE DEXCOM G6 CONTINUOUS GLUCOSE MONITORING SYSTEM



**“NO FINGERPRICKS?
BEST. NEWS. EVER!”**

Dexcom G6 eliminates fingerpricks for calibration and diabetes treatment decisions.*

* If your glucose alerts and readings from Dexcom G6 do not match symptoms or expectations, use a blood glucose meter to make diabetes treatment decisions.



MEET DEXCOM G6

Dexcom G6 provides real-time glucose readings every five minutes for people living with type 1 diabetes. Dexcom G6 lets you know where you've been, where you are and where you're headed!



Zero Fingerpricks & No Scanning Required*

Factory calibrated to allow for zero fingerpricks and no calibrations to make treatment decisions!

* If your glucose alerts and readings from the G6 do not match symptoms or expectations, use a blood glucose meter to make diabetes treatment decisions.



Longer-life Sensor†

Indicated to wear for up to 10 days



Slimmer‡, Water-resistant‡ Device

A slim, water-resistant sensor that is discreet and easy to wear



Urgent Low Soon Alert

Provides a 20-minute advance warning ahead of an urgent low



Customisable Alerts & Alarms

Set your range and get notified when you're heading high or low with 2 groups of alert schedules



Outstanding Accuracy

Even greater accuracy than our previous generation, offering you reliable results¹

† Compared to the Dexcom G5 Mobile CGM System.

‡ The Dexcom G6 Sensor and Transmitter are water resistant and may be submerged under 2.4 metres of water for up to 24 hours without failure when properly installed.




DEXCOM G6 FEATURES NOW WITH MORE PEACE OF MIND

Advanced Alert Options

- Urgent Low Soon alert - can provide a 20-minute advance warning of an urgent low (3.1 mmol/L)
- Alert Schedule - customise two groups of alerts in a 24 hour period (e.g. day and night schedules) for maximum flexibility

Updated CGM System Components

- Auto-applicator - simple, one-touch sensor insertion
- Discreet, low-profile device - 28% slimmer than Dexcom G5 Mobile CGM, making it nearly unnoticeable under clothing
- 10-day sensor wear - 43% longer wear than Dexcom G5 Mobile



REAL-TIME GLUCOSE READINGS

Dexcom G6 gives you the flexibility and convenience to monitor your glucose data when, where and how you want.

HOW DOES DEXCOM G6 WORK?



Applicator

- Simple, auto-applicator to insert the sensor just beneath the skin



Sensor

- Measures glucose levels for up to ten days
- Greater accuracy due to no paracetamol interference*



Transmitter

- Affixed to the top of sensor; sends data wirelessly to a display device



Display Device

Compatible Smart Device†

- Whether at work, home, or school you can check your glucose values right on your smart device with the Dexcom G6 app

Dexcom Receiver†

- Touchscreen display
- Displays a glucose reading every 5 minutes – up to 288 a day
- Trend arrows and vivid colours to show the speed and direction your glucose is heading

* Dexcom G6 readings can be used to make diabetes treatment decisions when taking up to a maximum paracetamol dose of 1,000mg every 6 hours. Taking a higher dose may affect readings.


† Smart device and Receiver not included with Dexcom G6 CGM system. For a list of compatible devices, please visit dexcom.com/compatibility

A CLOSER LOOK AT THE DEXCOM G6 MOBILE APP

The Dexcom G6 app displays glucose information at a glance. Customisable alerts and alarms can notify you of impending highs and lows.



The Dexcom G6 app is compatible with both iOS and Android devices.*



Want to view what the Dexcom G6 looks like?

Download the Dexcom G6 Simulator app.

This app displays samples of CGM scenarios, trends and various setup features associated with the Dexcom G6 CGM App.^




* For a list of compatible devices, please visit dexcom.com/compatibility

^ This app is for simulation purposes only and contains no live CGM data

THE POWER OF SHARING

With the built-in Share feature, up to 5 people can follow your glucose levels, giving you the circle of support you need.*

By downloading the Dexcom Follow App†, followers can view your glucose data directly from their compatible smart device‡, whether they live down the street or across the country.



* Do not make treatment decisions based upon Share and Follow readings. Always confirm with your compatible smart device or Receiver.

† Internet connectivity required to access Dexcom Follow.

‡ For a list of compatible devices, please visit dexcom.com/compatibility

MANAGE YOUR DIABETES WITH DEXCOM CLARITY

Dexcom
CLARITY

Dexcom CLARITY is an important part of your Dexcom CGM system. Using your CGM data, Dexcom CLARITY highlights glucose patterns, trends and statistics. You can share data with your clinic and monitor improvements between visits. Get to know CLARITY after you set up your Dexcom CGM.

WHEN AT HOME



Log in at clarity.dexcom.eu


Use your existing Dexcom login (if you have one).

WHILE ON THE GO



Get weekly notifications with the Dexcom CLARITY app

Notifications available to Dexcom mobile users.



NEXT STEPS

Follow the steps below to get started on Dexcom G6!

1

Learn more about Dexcom G6

Visit amsldiabetes.com.au to learn more about Dexcom. Speak to your healthcare professional to enquire if Dexcom G6 is right for you

2

Talk to AMSL Diabetes

For more information or any questions you may have, contact our 24/7 Australian-based customer care team at [1300 851 056](tel:1300851056) or email diabetes@amsl.com.au

3

Ready to order?

Learn about our Dexcom All-Rounder Subscription Program or our limited time Dexcom Starter Kit Offer. Once you're ready to place an order, visit us online at amsldiabetesshop.com.au to purchase


4

Start your Dexcom!

Once your box arrives, open your Dexcom G6 CGM system and follow the instructions. If you have any questions, please call our customer care team on [1300 851 056](tel:1300851056)

WHY AMSL DIABETES?


In addition to our 24/7 Australian-based expert customer care team, you get the benefit of:



Local Support



Shopping Online



Expert Training



Online Resources



Staying Updated



For any technical assistance or for information on any of our products, contact our Customer Care Team available on 1300 851 056.

For more information on Dexcom G6, please visit our website or contact us on **1300 851 056** or at **diabetes@amsl.com.au**

amsldiabetes.com.au    

Dexcom is only indicated for use in patients with type 1 diabetes over 2 years old and is not indicated for use in pregnancy or patients on dialysis treatment. References: 1. Wadwa RP, Laffel LM, Shah VN, Garg SK. Accuracy of a Factory-Calibrated, Real-Time Continuous Glucose Monitoring System During 10 Days of Use in Youth and Adults with Diabetes. Diabetes Technol Ther. 2018. Dexcom, Dexcom G6, Dexcom G5 Mobile, Dexcom G4 PLATINUM, Dexcom Share, Dexcom Follow and Dexcom CLARITY are registered trademarks of Dexcom, Inc. in the United States and/or other countries. ARTG 330535, PR-100-367 August 2021

