







# DEXCOM G6 & DEXCOM G5 MOBILE

## CONTINUOUS GLUCOSE MONITORING (CGM) SYSTEMS

	Dexcom G5 Mobile	Dexcom G6
<b>Calibration with fingerprick</b>	Twice in the first 2 hours then once every 12 hours	No calibration needed (0 fingerpricks required)*
<b>Sensor (including applicator)</b>		 1-button simplified insertion
<b>Sensor duration</b>	7 days	10 days
<b>Transmitter</b> Both G5 Mobile and G6 transmitters last 3 months		 28% lower profile than G5 Mobile
<b>Display option</b>	Smartphone (iOS, Android), smart watch†	Smartphone (iOS, Android), smart watch†
<b>Receiver (optional)</b>		 Intuitive touchscreen receiver
<b>Customisable real-time low/high BG alerts</b>	Yes	Yes • New "Urgent Low Soon" alert • New customisable day and nighttime alert schedules
<b>Paracetamol interference</b>	Yes	No‡
<b>Remote monitoring</b> §	Yes, up to 5 followers	Yes, up to 5 followers

\* If your glucose alerts and readings from the G6 do not match symptoms or expectations, use a blood glucose meter to make diabetes treatment decisions.

† To view a list of compatible devices, visit [www.dexcom.com/dexcom-international-compatibility](http://www.dexcom.com/dexcom-international-compatibility)

‡ If your glucose alerts and readings do not match your symptoms or expectations or you are taking medications containing paracetamol, you should perform a fingerprick to confirm your blood glucose level. A fingerprick is needed every 12 hours to calibrate.

§ Dexcom G6 readings can be used to make diabetes treatment decisions when taking a maximum paracetamol dose of 1000mg every 6 hours. Taking a higher dose may affect the readings.

¶ Requires the Dexcom Follow App and an internet connection. Followers should always confirm readings on the Dexcom G5 Mobile App, Dexcom G6 App or Receiver before making treatment decisions.

For more information on Dexcom G6, please contact us on **1300 851 056**  
or at [diabetes@amsl.com.au](mailto:diabetes@amsl.com.au)

[amsldiabetes.com.au](http://amsldiabetes.com.au)



© 2020 Dexcom, Inc. All rights reserved. Dexcom, Dexcom G6, Dexcom G5 Mobile are registered trademarks of Dexcom, Inc. in the United States and/or other countries. ARTG 330535, PR-100-378 April 2020

**DEXCOM**

