

t:slim X2™ Insulin Pump

T:LOCK™ INFUSION SETS AND CARTRIDGE

AutoSoft 90

- 90° insertion angle
- Integrated inserter
- Teflon cannula




AutoSoft 30

- 30° insertion angle
- Integrated inserter
- Teflon cannula



TruSteel Set

- Manual 90° insertion angle
- Stainless steel/metal cannula



VariSoft Set

- Manual angled set
- Teflon cannula



t:lock Infusion Set & Cartridge NDSS Order Codes (10 per pack)

AutoSoft 90 Infusion Set t:lock 6mm/60cm	826	TruSteel Infusion Set t:lock 6mm/60cm	834
AutoSoft 90 Infusion Set t:lock 6mm/110cm	827	TruSteel Infusion Set t:lock 6mm/80cm	835
AutoSoft 90 Infusion Set t:lock 9mm/60cm	828	TruSteel Infusion Set t:lock 8mm/60cm	836
AutoSoft 90 Infusion Set t:lock 9mm/110cm	829	TruSteel Infusion Set t:lock 8mm/80cm	837
AutoSoft 30 Infusion Set t:lock 13mm/60cm	830		
AutoSoft 30 Infusion Set t:lock 13mm/110cm	831		
VariSoft Infusion Set t:lock 13mm/60cm	832	t:lock Cartridge	
VariSoft Infusion Set t:lock 17mm/60cm	833	t:lock Cartridge	722

Cannula Fill and Site Change Guidelines for Infusion Sets

	Fill	Site Change
AutoSoft 90	6mm: 0.3 units, 9mm: 0.5 units	2-3 days
AutoSoft 30	0.7 units	2-3 days
TruSteel	0 units/no fill	1-2 days
VariSoft	0.7 units	2-3 days

1300 851 056
diabetes@amsl.com.au

amsl.com.au



TANDEM
Diabetes Care



a Dexcom company

ALWAYS READ THE LABEL AND FOLLOW THE DIRECTIONS FOR USE. Read the warnings available on amsl.com.au/resources before purchasing. Consult your healthcare professional to see which product is right for you.